

**Composition:**

Active ingredient: Hydrogen peroxide 3.0 %.

Inactive ingredients: Silver nitrate, fragrance, Aqua.

**Direction for use:**

Use one wipe to disinfect the surface. On heavily soiled surfaces, more than one wipe may be required. Allow to dry naturally. Do not use on patients or intensive medical devices. After use allow the surface to dry naturally before your next procedure. Do not re-use wipes.

**Direction for dispensing:**

- 1) Remove lid by lifting upward.
- 2) Pull towel from centre of the roll and thread through lid hole on the underside of lid.
- 3) Close lid. Pull towel from canister at 90 ° angle.
- 4) Snap hinged cap back into place to ensure towels stay moist.

**Storage:** Keep pack tightly closed when not in use. Store in cool, dry & ventilated place.

**Disposal:** Offer empty containers for recycling. If recycling is not available, discard container in trash. Discard towel in trash after use. Do not flush.

Note: For M.R.P., Batch No., Manufacturing date and Expiry date please see the lid.



**60 Rolled Wet Wipes**  
**6" x 8" Each**

# L-Hygiene™

## Hydrogen Peroxide 3% Wipes

### Silver Accelerated Hydrogen Peroxide

**Lint-free & tear resistant**  
**WHPNF00060**

**Ideal for disinfecting surfaces**

Home, Food service areas, Day care centers, Schools, Health care facilities, Universities and research laboratories, Wipe downs on filling and packing machinery, Restaurants and more

**Hard Surface Compatibility**

Baked enamel, Glass, Glazed ceramic, Glazed porcelain, Glazed tile, Laminated surfaces, Stainless Steel, and more

**Features & Benefits:**

Non- woven Polycellulose material

Each 6" x 8" (150 mm x 200 mm)

60 wipes/roll – MORE WIPES AND LESS WASTAGE

- Effective against common pathogens and more.
- No-Toxic Fumes
- Strongest disinfecting
- High Absorbency
- Plastic-safe

**Not Recommended For:**

Aluminium, Copper, Galvanized steel, Silver, Rubber, Porous plastic, Wood

**NON STERILE-NOT INTENDED FOR MEDICAL USE****KEEP OUT OF REACH OF CHILDREN**

Do not use in or near the eyes. In case of contact, rinse eyes thoroughly with water. Discontinue use if irritation develops.



**WARNINGS:** For external use only.

Manufactured by:

**Loba Chemie Pvt. Ltd.**

Add.: Plot No.- D 22, M.I.D.C., Tarapur, Boisar,  
Palghar - 401506, Maharashtra(India)

Customer Care No.: +91 22 66636663

Email: [info@l-hygiene.com](mailto:info@l-hygiene.com)

Web: [www.l-hygiene.com](http://www.l-hygiene.com)